



BIO•PHARMA S E R V I C E S

— THE POWER TO MAKE® —

CORYNEX PROTEIN EXPRESSION:
SIMPLIFYING RECOMBINANT PROTEIN PRODUCTION



The rise of biologics is having a dramatic impact on how physicians treat many serious and chronic illnesses. While these drugs have great promise, biologics are complex and, as a result, are very expensive to manufacture and subject to technical pitfalls. Pharmaceutical companies looking to add biologics to their pipelines must come up with innovative ways to approach drug development and cut the high costs associated with it.

One way to do this is to simplify the processes to produce and purify recombinant proteins. Simpler processes result in reduced production costs and a faster time to market. In this white paper, we review how the Corynex[®] Protein Expression System uses an extensive toolbox to improve the levels of protein secretion and overall success rate.

THE POWER TO MAKE[®]

WHAT ARE SOME OF THE BIGGEST CHALLENGES OF PRODUCING AND PURIFYING RECOMBINANT PROTEINS?

Manufacturing highly purified recombinant proteins requires 1) an upstream with high titer and recovery of active, correctly folded proteins, 2) a downstream with multiple steps of purification with high yield, and, therefore, 3) significant time and cost. Recombinant proteins are often expressed in *E. coli* as insoluble aggregates called inclusion bodies. Active proteins can be obtained after dissolving the inclusion bodies and refolding the protein, which generally returns quite low yield. Expression using mammalian systems requires long cycles for cell breeding, screening, and cultivation. Further, even a single amino acid substitution can cause significantly different levels of expression in any system. This results in researchers using trial and error approaches to look for the appropriate expression system for each of their proteins of interest.

WHAT IS CORYNEX?

Developed by Ajinomoto Co., Inc., Corynex® is an innovative protein/peptide expression system using a gram-positive, non-sporulating soil bacterium, *Corynebacterium glutamicum*. This nonpathogenic bacterium has successfully been applied in the industrial production of amino acids utilized in human food, animal feed, and pharmaceutical products for more than 60 years.

Due to this rich history, the metabolic physiology, biochemistry, genetics, and fermentation systems of *C. glutamicum* are well known, and various technologies have been developed in order to improve the productivity of amino acids produced with this system. Based on this knowledge, *C. glutamicum* has been developed not only as an

amino acid producing host, but also as a protein production host named Corynex.

Corynex has the advantage of a scalable, high-cell-density fermentation process that can be used to manufacture commercially valuable proteins such as biopharmaceuticals, drug targets, and enzymes. Corynex secretes soluble, properly-folded, and biologically-active recombinant proteins directly into the culture media with high purity.

WHAT MAKES CORYNEX DIFFERENT FROM OTHER SIMILAR TECHNOLOGIES?

Microbial expression systems, such as *E. coli*, are inexpensive and easy to handle in upstream, although in most cases, require refolding and removal of host cell proteins and endotoxins in complex downstream operations. On the other hand, mammalian cells, such as CHO cells, have highly sophisticated protein modification and maintenance mechanisms, which make downstream processes simpler than many microbial systems. However, creation of stable cell lines and bioreactor production require substantial time and cost. Corynex possesses the best of both worlds: **simplified upstream and downstream processes.**

HOW IS THE PURIFICATION PROCESS USED BY CORYNEX DIFFERENT FROM TRADITIONAL TECHNOLOGIES?

Corynex has the ability to:

- **Secrete with high purity** (Figure 1): Because target proteins/peptides are directly secreted into the culture media, there is no need for cell disruption. Purity of the secreted protein is quite high since *Corynebacterium* secretes a limited amount of host cell proteins.

Furthermore, *C. glutamicum* produces no endotoxins. These characteristics make it possible to construct simple processes with reduced purification steps. High purity of secreted proteins/peptides can be observed in SDS-PAGE analyses of culture supernatant from expression tests.




FIGURE 1: High purity of secreted proteins (unpurified culture supernatants)

- **No extracellular protease degradation:** Protease activities are not detected in the culture supernatant, and the secreted protein is not degraded, which is often a problem in many other protein secretion systems. Corynex has succeeded in secretion expression of various proteins, such as MTG (microbial transglutaminase), epidermal growth factor (EGF) derived from human, insulin-like growth factor (IGF -1), basic fibroblast growth factor (bFGF), and VHH antibody, however, degradation by proteases derived from the host was not observed (Figure 1). This property is also an advantage in the construction of simple purification processes.
- **Express active proteins (Figure 2):** Proteins having complex disulfide bonds (S-S bonds) can be secreted in active form using the Corynex system. There is no need for inclusion system. There is no need for inclusion body recovery, or

complicated and unwanted refolding processes. In the expression of human epidermal growth factor (hEGF), which has three intramolecular S-S bonds, it was secreted as a single band in SDS-PAGE analysis, and confirmed to be properly folded and functional with correct S-S bonds. The hEGF produced using the Corynex system had similar biological activity as endogenous hEGF. In addition, not only were proteins with intramolecular S-S bonds successfully secreted into active form, but also the heterodimer antibody fragment Fab with intermolecular S-S bonds.

Because there are many fewer steps than required with other bacterial platforms, there is less protein loss during purification and overall yield is increased. This results in a reduction of costs and time in the development of protein therapeutics.

WHAT SYSTEMS ARE USED IN CORYNEX?

As with all other recombinant protein production platforms, Corynex cannot universally secrete all proteins. To address this issue, three secretion technologies were developed as a standard system in initial expression trials.

- **Sec pathway:** The first technology is the general secretory pathway called “Sec system”. The Sec system is a widely conserved pathway in both prokaryotes and eukaryotes, in which translated proteins are secreted through the Sec machinery in unfolded form, and are subsequently folded extracellularly.
- **CspB fusion method:** The second is “CspB fusion method” which is an application technology of Sec-type secretion. CspB is a cell surface protein secreted in large quantities in *C. glutamicum*. Some proteins can be easily

secreted by fusing the N-terminal sequence residues of CspB with the N-terminal of the target protein as a tag. The intact protein can be obtained by using a specific protease which recognizes the sequence design between CspB-tag and target protein. This is a useful method when the sufficient secretion cannot be obtained by the basic Sec system.

- **Tat pathway:** The third technology is “Tat system”, a pathway quite distinct from the Sec system. The more recently identified Tat secretory pathway is able to transport folded proteins, including those that are large and heterologous, which are sometimes poorly secreted by the Sec pathway (Figure 2). The Tat system is generally considered to be a route for secretion of small amounts of proteins, but in Corynex, it has been successfully engineered to secrete various recombinant proteins with high productivity.

Secretion of many properly folded recombinant proteins has been achieved with these pathways including hundreds with disulfide bonded and dimerized structures.




FIGURE 2: Illustration of the Corynex system’s three expression constructs.

WHAT TOOLS ARE USED WITH CORYNEX?

Since it is not always possible to obtain high titers in the initial expression trial, a toolbox has been developed for Corynex to achieve higher

productivity during strain and process development (Figure 3).

- **Signal peptide library:** A library of signal peptides derived from *C. glutamicum* can be screened to further increase secretion levels. The optimum signal peptide is determined empirically for each target protein. Secretion levels of the recombinant protein could be increased by screening for optimal signal peptides.
- **High performance host library:** A library of mutant hosts with specific gene deficiencies or mutations by random mutagenesis, which often improve productivity. It is possible to find a better host in comparison with the standard host for each target protein.
- **Tailored mutagenesis:** Nitrosoguanidine (NG)-induced mutagenesis and high throughput screening allows the Corynex team to quickly and efficiently select optimal host strains tailored to each target protein. More than 2,000 strains can be generated and tested within three months. With optimized host strains, protein expression for target proteins can be significantly increased.




FIGURE 3: Toolboxes for strain development.

FERMENTATION PROCESS DEVELOPMENT

The fermentation process using Corynex is quite simple. As previously mentioned, *C. glutamicum* has been used in the industrial production of amino acids for 60 years. Therefore, the cultivation method using a simple synthetic medium has been well-established, and there is a lot of know-how for process development.

As shown in Figure 4, the combination of strain and fermentation process development can drastically improve the productivity of VHH antibodies.




FIGURE 4: Example of VHH antibody production improvement.

WHAT IS THE TRACK RECORD OF THE CORYNEX PLATFORM?

Biologics candidates developed by global pharmaceutical companies using Corynex have entered Phase I clinical trials. After feasibility studies were successfully performed, the pharma company engaged with our team to scale-up and produce GMP material at the 1000-liter scale to support Phase 1 clinical studies. Successful GMP batches have been produced and the first human subjects have been dosed.

As an industry leader in microbial fermentation, Ajinomoto Bio-Pharma Services' expertise and knowledge of incorporating the proprietary Corynex system to simplify the purification process, results in improved delivery of complex recombinant proteins for our clients. By partnering with a CDMO with extensive large molecule experience who can utilize innovative platforms, such as Corynex, you will find yourself with an improved drug development process, reduced production costs and faster time to market.



info@US.AjiBio-Pharma.com

www.AjiBio-Pharma.com